


Anti-Dhtkd1 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Hy7

Item Number: CM0025527



Introduction

Rabbit polyclonal antibody targeting human DHTKD1 (2-oxoadipate dehydrogenase E1 component), validated for WB and ELISA. Cross-reacts with human, mouse, rat; detects ~103 kDa band. Ideal for mitochondrial metabolism and lysine/tryptophan catabolism research.

[Learn More](#)

Parameter	Value
Product Name	DHTKD1 Rabbit pAb
Remarks/Alias	E1a; AAKAD; CMT2Q; AMOXAD; OADC-E1; OADH-E1; DHTKD1
Species	Human
GeneID (Human)	55526
GeneID	55526
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human DHTKD1 (NP_061176.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96HY7
Protein Weight	103kDa
Shipping	Ice bag