


Anti-Dhx30 Rabbit Polyclonal Antibody For Wb, Elisa - Q7L2E3

Item Number: CM0026110



Introduction

Rabbit polyclonal antibody targeting DHX30, suitable for WB and ELISA. Recognizes human, mouse, and rat DHX30. Supports research on mitochondrial ribosome assembly and neurogenesis.

[Learn More](#)

Parameter	Value
Product Name	DHX30 Rabbit pAb
Remarks/Alias	DDX30; RETCOR; NEDMIAL; DHX30
Species	Human
GeneID (Human)	22907
GeneID	22907
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 906-1155 of human DHX30 (NP_055781.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7L2E3
Protein Weight	134kDa
Shipping	Ice bag