


Anti-Dhx30 Rabbit Polyclonal Antibody For Wb, Elisa - Q7L2E3

Item Number: CM0022660



Introduction

Rabbit polyclonal anti-DHX30 antibody for WB and ELISA, detecting human and mouse DHX30 (~134kDa). DHX30 is a DEAH-box RNA helicase involved in mitochondrial ribosome biogenesis, neurogenesis, and innate immunity, making this antibody suitable for mitochondrial and neuroscience research.

[Learn More](#)

Parameter	Details
Product Name	DHX30 Rabbit pAb
Remarks/Alias	DDX30, RETCOR, NEDMIAL, DHX30
Species	Human
Gene ID (Human)	22907
Gene ID	22907
Immunogen	Recombinant protein corresponding to a sequence within amino acids 336-501 of human DHX30 (NP_619520.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7L2E3
Protein Weight	134kDa
Shipping	Ice bag