


Anti-Dhx36 Polyclonal Antibody For Wb, Ip, Elisa - Q9H2U1

Item Number: CM0030296



Introduction

DHX36 Rabbit Polyclonal Antibody validated for WB, IP, and ELISA. Reacts with human, mouse, and rat DHX36, a G-quadruplex resolving helicase involved in gene regulation and genomic integrity.

[Learn More](#)

Parameter	Value
Product Name	DHX36 Rabbit pAb
Remarks/Alias	G4R1; RHAU; DDX36; MLEL1; DHX36
Species	Human
Gene ID (Human)	170506
Gene ID	170506
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human DHX36 (NP_065916.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H2U1
Protein Weight	115kDa
Shipping	Ice bag