

Anti-Dhx38 Rabbit Polyclonal Antibody For Wb, Elisa - Q92620

Item Number: CM0025191

Anti-DHX38 Rabbit Polyclonal Antibody for WB, ELISA - Q92620

Introduction

DHX38 Rabbit Polyclonal Antibody validated for Western Blot and ELISA. Reacts with human, mouse, and rat. Targets pre-mRNA splicing factor DHX38, a spliceosome component and ATP-dependent RNA helicase.

[Learn More](#)

Parameter	Value
Product Name	DHX38 Rabbit pAb
Remarks/Alias	RP84, DDX38, PRP16, PRPF16, DHX38
Species	Human
GeneID (Human)	9785
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1133-1227 of human DHX38 (NP_054722.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92620
Protein Weight	141kDa
Shipping	Ice bag