

Anti-Dhx57 Polyclonal Antibody For Wb And Elisa - Q6P158

Item Number: CM0021868

Introduction

Rabbit polyclonal antibody targeting human DHX57/DDX57, suitable for Western blot (WB) and ELISA. Cross-reacts with mouse. For studies on ATP-dependent RNA helicase function.

[Learn More](#)

Parameter	Value
Product Name	DHX57 Rabbit pAb
Remarks/Alias	DDX57; DHX57
Species	Human
GeneID (Human)	90957
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 700-900 of human DHX57 (NP_945314.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q6P158
Protein Weight	156kDa
Shipping	Ice bag