

Anti-Dhx58 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96C10

Item Number: CM0028282



Introduction

Rabbit polyclonal antibody against DHX58 (LGP2), validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. For innate immunity and antiviral signaling research.

[Learn More](#)

Parameter	Value
Product Name	DHX58 Rabbit pAb
Remarks/Alias	LGP2, RLR-3, D11LGP2, D11lgp2e, DHX58
Species	Human
GeneID (Human)	79132
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 389-678 of human DHX58 (NP_077024.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96C10
Protein Weight	77kDa
Shipping	Ice bag