

Anti-Dhx8 Rabbit Polyclonal Antibody For Wb, Elisa - Q14562

Item Number: CM0015839



Introduction

Human DHX8 rabbit polyclonal antibody validated for Western blot and ELISA. Targets the ATP-dependent RNA helicase involved in pre-mRNA splicing and nuclear export. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	DHX8 Rabbit pAb
Remarks/Alias	DDX8, Dhr2, HRH1, PRP22, PRPF22, DHX8
Species	Human
GeneID (Human)	1659
GeneID	1659
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1001-1220 of human DHX8 (NP_004932.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14562
Protein Weight	139kDa
Shipping	Ice bag