


Anti-Dhx9/Rna Helicase A Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q08211

Item Number: CM0020636



Introduction

Rabbit polyclonal antibody targeting DHX9/RNA Helicase A, suitable for WB, IF/ICC, IP, and ELISA. Reacts with human, mouse, and rat samples. Ideal for research on transcriptional regulation, mRNA processing, and viral infection.

[Learn More](#)

Parameter	Specification
Product Name	DHX9/RNA Helicase A Rabbit pAb
Remarks/Alias	LKP; RHA; DDX9; NDH2; NDH1; DHX9/RNA Helicase A
Species	Human
GeneID (Human)	1660
GeneID	1660
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08211
Protein Weight	141kDa
Shipping	Ice bag