


Anti-Dhx9/Rna Helicase A Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q08211

Item Number: CM0012487



Introduction

Rabbit monoclonal anti-DHX9/RNA Helicase A antibody validated for WB, IHC-P, IF/ICC, IP, ELISA. Recognizes human, mouse, rat DHX9 (141 kDa). Ideal for transcription, RNA processing, and helicase studies.

[Learn More](#)

Parameter	Value
Product Name	DHX9/RNA Helicase A Rabbit mAb
Remarks/Alias	LKP; RHA; DDX9; NDH2; NDH3; DHX9/RNA Helicase A
Species	Human
GeneID (Human)	1660
GeneID	1660
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1100-1200 of human DHX9/RNA Helicase A (Q08211)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08211
Protein Weight	141kDa
Shipping	Ice bag