


Anti-Diaph1 Polyclonal Antibody (Rabbit) For Wb And Elisa - O60610

Item Number: CM0022693



Introduction

Rabbit anti-DIAPH1 polyclonal antibody validated for Western blot (WB) and ELISA. Targets human DIAPH1 (141 kDa), an actin nucleation factor involved in cytoskeleton organization. Ideal for cell biology and hearing research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] DIAPH1 Rabbit pAb
Remarks/Alias	DIA1; DRF1; DFNA1; LFHL1; SCBMS; hDIA1; mDia1; [KD Validated] DIAPH1
Species	Human
GeneID (Human)	1729
GeneID	1729
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1158-1272 of human DIAPH1 (NP_005210.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O60610
Protein Weight	141kDa
Shipping	Ice bag