


Anti-Diaph1 Rabbit Monoclonal Antibody For Wb And Elisa - O60610

Item Number: CM0002907



Introduction

Rabbit monoclonal antibody specific for human DIAPH1 (DIA1/DRF1). Validated for WB and ELISA; cross-reacts with mouse and rat. Suitable for actin cytoskeleton and cell migration studies.

[Learn More](#)

Parameter	Value
Product Name	DIAPH1 Rabbit mAb
Remarks/Alias	DIA1; DRF1; DFNA1; LFHL1; SCBMS; hDIA1; mDia1; DIAPH1
Species	Human
GeneID (Human)	1729
GeneID	1729
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human DIAPH1 (O60610)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60610
Protein Weight	141 kDa
Shipping	Ice bag