

Anti-Diaph2 Rabbit Polyclonal Antibody For Wb, Elisa - O60879

Item Number: CM0017361

Introduction

DIAPH2 rabbit polyclonal antibody validated for WB and ELISA. Detects DIAPH2 in Human, Mouse, Rat samples. Relevant for studying oogenesis, endosomal trafficking, and actin cytoskeleton regulation.

[Learn More](#)

Parameter	Specification
Product Name	DIAPH2 Rabbit pAb
Remarks/Alias	DIA; POF; DIA2; DRF2; POF2; POF2A; DIAPH2
Species	Human
Gene ID (Human)	1730
Gene ID	1730
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human DIAPH2 (NP_009293.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60879
Protein Weight	126kDa
Shipping	Ice bag