


Anti-Diaph3 Polyclonal Antibody For Wb, Elisa - Q9Nsv4

Item Number: CM0016574



Introduction

Human DIAPH3 Rabbit Polyclonal Antibody validated for WB and ELISA. Targets diaphanous-related formin-3 (mDia2) involved in actin cytoskeleton dynamics, cytokinesis, and SRF activation. Suitable for research on hearing loss.

[Learn More](#)

Parameter	Value
Product Name	DIAPH3 Rabbit pAb
Remarks/Alias	AN; DIA2; DRF3; AUNA1; NSDAN; diap3; mDia2; DIAPH3
Species	Human
GeneID (human)	81624
GeneID	81624
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NSV4
Protein Weight	137kDa
Shipping	Ice bag