

Anti-Dicer1 Polyclonal Antibody For If/Icc, Elisa - Q9Upy3

Item Number: CM0020129



Introduction

DICER1 rabbit polyclonal antibody for human samples, validated for IF/ICC and ELISA. Targets endoribonuclease Dicer, key in miRNA and siRNA processing.

[Learn More](#)

Parameter	Description
Product Name	DICER1 Rabbit pAb
Remarks/Alias	DCR1; GLOW; MNG1; aviD; Dicer; HERNA; RMSE2; Dicer1e; K12H4.8-LIKE; DICER1
Species	Human
Gene ID (Human)	23405
Gene ID	23405
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1100-1200 of human DICER1 (NP_001258211.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9UPY3
Protein Weight	219kDa
Shipping	Ice bag