

Anti-Digoxigenin/Digoxin Rabbit Monoclonal Antibody For Db And Elisa Applications

Item Number: CM0000004

Introduction

Rabbit monoclonal antibody targeting digoxigenin/digoxin. Suitable for DB and ELISA with species-independent reactivity. Ideal for IVD assay development and digoxin monitoring.

[Learn More](#)

Parameter	Value
Product Name	Anti-Digoxigenin/Digoxin Rabbit mAb
Immunogen	Chemical compounds Digoxin
Source	Rabbit
Category	Other Modified Antibodies
Application	DB, ELISA
Cross Reactivity	Species independent
Shipping	Ice bag