


Anti-Dimethyl-Histone H3 (Lys79) Rabbit Mab For Chromatin/Epigenetic Research (Q16695/P68431)

Item Number: CM0000354



Introduction

Rabbit monoclonal antibody specific for dimethylated histone H3 at lysine 79 (H3K79me2). Validated for WB, ChIP-seq, IF, ELISA in human, mouse, rat. Ideal for epigenetic and chromatin analysis.

[Learn More](#)

Parameter	Value
Product Name	DiMethyl-Histone H3-K79 Rabbit mAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; DiMethyl-Histone H3-K79
Species	Human
GeneID (Human)	8290/ 8350
GeneID	8290/8350
Immunogen	Synthetic peptide
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, DB, ChIP, ChIP-seq, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag