


Anti-Dimt1 Polyclonal Antibody For Wb, Elisa - Q9Unq2

Item Number: CM0030607



Introduction

High-quality rabbit polyclonal antibody targeting DIMT1, a key enzyme in 18S rRNA dimethylation and ribosome biogenesis. Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Ideal for nucleolar and rRNA processing studies.

[Learn More](#)

Parameter	Value
Product Name	DIMT1 Rabbit pAb
Remarks/Alias	DIM1; DIMT1L; HUSSY5; HSA9761; DIMT1
Species	Human
GeneID (Human)	27292
GeneID	27292
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-126 of human DIMT1 (NP_055288.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UNQ2
Protein Weight	35kDa
Shipping	Ice bag