

Anti-Dio3 Rabbit Polyclonal Antibody For Wb & Elisa - P55073

Item Number: CM0027174

Introduction

Rabbit polyclonal antibody targeting human DIO3 (Thyroxine 5-deiodinase), suitable for WB and ELISA applications. Ideal for studying thyroid hormone metabolism and related research - P55073.

[Learn More](#)

Parameter	Value
Product Name	DIO3 Rabbit pAb
Remarks/Alias	D3; 5DIII; TXDI3; DIOIII; DIO3
Species	Human
GeneID (Human)	1735
GeneID	1735
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 70-170 of human DIO3 (NP_001353.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P55073
Protein Weight	34kDa
Shipping	Ice bag