

Anti-Dip2C Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y2E4

Item Number: CM0022704

Introduction

DIP2C rabbit polyclonal antibody validated for WB, IF/ICC, ELISA with reactivity against human, mouse, and rat. Recombinant immunogen corresponding to amino acids 100-260.

[Learn More](#)

Parameter	Value
Product Name	DIP2C Rabbit pAb
Remarks/Alias	KIAA0934; DIP2C
Species	Human
GeneID (Human)	22982
GeneID	22982
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-260 of human DIP2C (NP_055789.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2E4
Protein Weight	171kDa
Shipping	Ice bag