


Anti-Diras3 Rabbit Polyclonal Antibody For Wb, Elisa - 095661

Item Number: CM0031042



Introduction

DIRAS3 Rabbit Polyclonal Antibody validated for WB and ELISA. Recognizes human DIRAS3, cross-reacts with mouse and rat. Ideal for cancer research and tumor suppressor studies.

[Learn More](#)

Parameter	Value
Product Name	DIRAS3 Rabbit pAb
Remarks/Alias	ARHI; NOEY2; DIRAS3
Species	Human
GeneID (Human)	9077
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-229 of human DIRAS3 (NP_004666.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O95661
Protein Weight	26kDa
Shipping	Ice bag