

Anti-Disc1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nri5

Item Number: CM0028110



Introduction

High-quality rabbit polyclonal antibody against human DISC1 (Q9NRI5) validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for neuroscience and schizophrenia research.

[Learn More](#)

Parameter	Value
Product Name	DISC1 Rabbit pAb
Remarks/Alias	SCZD9; C1orf136; DISC1
Species	Human
Gene ID (Human)	27185
Gene ID	27185
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human DISC1 (NP_001158009.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRI5
Protein Weight	94kDa
Shipping	Ice bag