

Anti-Disp1 Rabbit Polyclonal Antibody For Wb, Elisa - Q96F81

Item Number: CM0017549



Introduction

Rabbit polyclonal antibody against human DISP1, validated for Western blot and ELISA. Targets protein dispatched homolog 1, a key regulator of hedgehog signaling. Ideal for developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	DISP1 Rabbit pAb
Remarks/Alias	DISPA; DISP1
Species	Human
GeneID (Human)	84976
GeneID	84976
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1325-1524 of human DISP1 (NP_116279.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96F81
Protein Weight	171kDa
Shipping	Ice bag