


Anti-DJ-1/Park7 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q99497

Item Number: CM0021458



Introduction

High-quality rabbit polyclonal antibody against human DJ-1/PARK7, validated for WB, IHC-P, IF/ICC, and ELISA, with cross-reactivity to mouse and rat. Ideal for oxidative stress, neurodegeneration, and deglycase research.

[Learn More](#)

Parameter	Value
Product Name	DJ-1/PARK7 Rabbit pAb
Remarks/Alias	DJ1; DJ-1; GATD2; HEL-S-67p; DJ-1/PARK7
Species	Human
Gene ID (Human)	11315
Gene ID	11315
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-189 of human DJ-1/PARK7 (NP_009193.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99497
Protein Weight	20kDa
Shipping	Ice bag