


## Anti-DJ-1/Park7 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q99497

Item Number: CM0003807



### Introduction

High-quality anti-DJ-1/PARK7 rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat DJ-1/PARK7 (20 kDa), a multifunctional oxidative stress protein. Ideal for neurodegeneration and cytoprotection research.

[Learn More](#)

| Parameter        | Value                                                                                                 |
|------------------|-------------------------------------------------------------------------------------------------------|
| Product Name     | DJ-1/PARK7 Rabbit mAb                                                                                 |
| Remarks/Alias    | DJ1; DJ-1; GATD2; HEL-S-67p; DJ-1/PARK7                                                               |
| Species          | Human                                                                                                 |
| GeneID (Human)   | 11315                                                                                                 |
| GeneID           | 11315                                                                                                 |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DJ-1/PARK7 (Q99497) |
| Source           | Rabbit                                                                                                |
| Category         | Monoclonal Antibodies                                                                                 |
| Application      | WB, IHC-P, ELISA                                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                                                     |
| SWISS            | Q99497                                                                                                |
| Protein Weight   | 20kDa                                                                                                 |
| Shipping         | Ice bag                                                                                               |