


Anti-Dkc1 Polyclonal Antibody For Ihc-P, If/Icc, Elisa - O60832

Item Number: CM0022719



Introduction

Rabbit polyclonal antibody against human DKC1 (dyskerin), validated for IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for ribosome biogenesis, telomere maintenance, and dyskeratosis congenita research.

[Learn More](#)

Parameter	Specification
Product Name	DKC1 Rabbit pAb
Remarks/Alias	DKC; CBF5; DKCX; NAP57; NOLA4; XAP101; DKC1
Species	Human
GeneID (human)	1736
GeneID	1736
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60832
Protein Weight	58kDa
Shipping	Ice bag