

Anti-Dkc1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O60832

Item Number: CM0008482



Introduction

High-quality anti-DKC1 rabbit monoclonal antibody validated for Western blot, IHC-P, and ELISA. Detects human, mouse, rat DKC1 (dyskerin), involved in pseudouridylation, ribosome biogenesis, and telomere maintenance.

[Learn More](#)

Parameter	Value
Product Name	DKC1 Rabbit mAb
Remarks/Alias	DKC; CBF5; DKCX; NAP57; NOLA4; XAP101; DKC1
Species	Human
Gene ID (Human)	1736
Gene ID	1736
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human DKC1 (O60832).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60832
Protein Weight	58kDa
Shipping	Ice bag