


Anti-Dkc1 Rabbit Polyclonal Antibody For Wb - 060832

Item Number: CM0017437



Introduction

Rabbit polyclonal antibody against human DKC1, for WB, IHC-P, IF/ICC, ELISA. Suitable for ribosome biogenesis and telomere maintenance research. Cross-reacts with human and mouse. Predicted molecular weight 58kDa.

[Learn More](#)

Parameter	Value
Product Name	DKC1 Rabbit pAb
Remarks/Alias	DKC; CBF5; DKCX; NAP57; NOLA4; XAP101; DKC1
Species	Human
GeneID (Human)	1736
GeneID	1736
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human DKC1 (NP_001354.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60832
Protein Weight	58kDa
Shipping	Ice bag