

Anti-Dkk1 Rabbit Polyclonal Antibody For Wb And Elisa - O94907

Item Number: CM0030960



Introduction

This high-quality anti-DKK1 rabbit polyclonal antibody is validated for WB and ELISA. It detects human, mouse, and rat DKK1, making it ideal for Wnt signaling and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	DKK1 Rabbit pAb
Remarks/Alias	SK; DKK-1; DKK1
Species	Human
Gene ID (Human)	22943
Gene ID	22943
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-266 of human DKK1 (NP_036374.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94907
Protein Weight	29kDa
Shipping	Ice bag