

Anti-Dkk3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ubp4

Item Number: CM0019812

Introduction

High-quality rabbit polyclonal antibody against human DKK3 (Dickkopf-related protein 3). Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for Wnt signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	DKK3 Rabbit pAb
Remarks/Alias	RIG; CRRL; REIC; DKK3
Species	Human
GeneID (Human)	27122
GeneID	27122
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-200 of human DKK3 (NP_056965.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBP4
Protein Weight	38kDa
Shipping	Ice bag