

Anti-Dkk3 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q9Ubp4

Item Number: CM0007689

Introduction

Recombinant anti-DKK3 rabbit monoclonal antibody validated for WB, IF-P, and ELISA. Detects human, mouse, and rat DKK3 protein, a Wnt signaling antagonist involved in development and disease.

[Learn More](#)

Parameter	Value
Product Name	DKK3 Rabbit mAb
Remarks/Alias	RIG; CRRL; REIC; DKK3
Species	Human
GeneID (Human)	27122
GeneID	27122
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-200 of human DKK3 (NP_056965.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBP4
Protein Weight	38kDa
Shipping	Ice bag