

Anti-Dkk3 Rabbit Monoclonal Antibody For Wb/If/Icc/Elisa - Q9Ubp4

Item Number: CM0007811

Introduction

Recombinant rabbit monoclonal antibody against human DKK3 (Q9UBP4), validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for Wnt signaling and developmental biology research.

[Learn More](#)

Parameter	Details
Product Name	DKK3 Rabbit mAb
Remarks/Alias	RIG; CRRL; REIC; DKK3
Species	Human
Gene ID (Human)	27122
Gene ID	27122
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-350 of human DKK3 (Q9UBP4) (NP_056965.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBP4
Protein Weight	38kDa
Shipping	Ice bag