

Anti-Dkk3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ubp4

Item Number: CM0022365

Introduction

Rabbit polyclonal antibody targeting human DKK3 (Dickkopf-related protein 3). Suitable for WB, IF/ICC, and ELISA applications. Cross-reacts with human, mouse, and rat. Ideal for Wnt signaling, development, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	DKK3 Rabbit pAb
Remarks/Alias	RIG; CRRL; REIC; DKK3
Species	Human
GenelD (Human)	27122
GenelD	27122
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-200 of human DKK3 (NP_056965.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBP4
Protein Weight	38kDa
Shipping	Ice bag