


Anti-Dkk4 (Dkk-4) Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ubt3

Item Number: CM0029736



Introduction

High-quality DKK4 rabbit polyclonal antibody validated for WB and ELISA. Cross-reacts with human and mouse DKK4. Ideal for studying Wnt signaling, development, and related diseases. Source your research materials with confidence.

[Learn More](#)

Parameter	Value
Product Name	DKK4 Rabbit pAb
Remarks/Alias	DKK-4; DKK4
Species	Human
Gene ID (Human)	27121
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-224 of human DKK4 (NP_055235.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UBT3
Protein Weight	25kDa
Shipping	Ice bag