


## Anti-Dlat Polyclonal Antibody For Wb, If/Icc, Elisa - P10515

Item Number: CM0028785



### Introduction

Rabbit polyclonal anti-DLAT antibody for Western blot, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Suitable for mitochondrial metabolism and pyruvate dehydrogenase complex research.

[Learn More](#)

| Parameter        | Value                              |
|------------------|------------------------------------|
| Product Name     | DLAT Rabbit pAb                    |
| Remarks/Alias    | E2; PBC; DLTA; PDCE2; PDC-E2; DLAT |
| Species          | Human                              |
| GeneID (human)   | 1737                               |
| GeneID           | 1737                               |
| Immunogen        | Recombinant protein                |
| Source           | Rabbit                             |
| Category         | Polyclonal Antibodies              |
| Application      | WB, IF/ICC, ELISA                  |
| Cross Reactivity | Human, Mouse, Rat                  |
| SWISS            | P10515                             |
| Protein Weight   | 69kDa                              |
| Shipping         | Ice bag                            |