


Anti-Dlat Rabbit Monoclonal Antibody For Wb, Elisa - P10515

Item Number: CM0009870



Introduction

Rabbit monoclonal antibody specific for DLAT (PDC-E2), the dihydrolipoamide acetyltransferase component of the pyruvate dehydrogenase complex. Validated for WB and ELISA, it cross-reacts with human, mouse, and rat. Ideal for mitochondrial metabolism research.

[Learn More](#)

Parameter	Value
Product Name	DLAT Rabbit mAb
Remarks/Alias	E2; PBC; DLTA; PDCE2; PDC-E2; DLAT
Species	Human
GeneID (human)	1737
GeneID	1737
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 548-647 of human DLAT (P10515)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10515
Protein Weight	69 kDa
Shipping	Ice bag