


Anti-Dlat Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P10515

Item Number: CM0010793



Introduction

Rabbit monoclonal antibody against human DLAT (dihydrolipoamide acetyltransferase) for multi-application use: WB, IF/ICC, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for mitochondrial metabolism and PBC research.

[Learn More](#)

Parameter	Value
Product Name	DLAT Rabbit mAb
Remarks/Alias	E2, PBC, DLTA, PDCE2, PDC-E2
Species	Human
GeneID (Human)	1737
GeneID	1737
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 87-270 of human DLAT (NP_001922.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10515
Protein Weight	69 kDa
Shipping	Ice bag