


Anti-Dlc1 Polyclonal Antibody For Wb And Elisa - Q96Qb1

Item Number: CM0022535



Introduction

CamelBio supplies rabbit anti-DLC1 polyclonal antibody (Q96QB1), validated for WB and ELISA. Detects human, mouse, rat DLC1. Recombinant immunogen (aa 204-463). Ideal for tumor suppressor and cell migration research.

[Learn More](#)

Parameter	Value
Product Name	DLC1 Rabbit pAb
Remarks/Alias	HP; ARHGAP7; STARD12; p122-RhoGAP; DLC1
Species	Human
GeneID (Human)	10395
GeneID	10395
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 204-463 of human DLC1 (NP_079043.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96QB1
Protein Weight	171kDa
Shipping	Ice bag