

Anti-Dld/Dldh Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P09622

Item Number: CM0016498



Introduction

Dihydrolipoamide dehydrogenase rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human and mouse DLD/DLDH. Suitable for studying mitochondrial metabolism, histone succinylation, and sperm capacitation.

[Learn More](#)

Parameter	Value
Product Name	DLDH/DLD Rabbit pAb
Remarks/Alias	E3; LAD; DLDD; DLDH; GCSL; PHE3; OGDC-E3; DLDH/DLD
Species	Human
GeneID (Human)	1738
GeneID	1738
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 36-300 of human DLDH/DLDH/DLDH/DLD (NP_000099.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09622
Protein Weight	54kDa
Shipping	Ice bag