

Anti-Dldh/Dld Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P09622

Item Number: CM0008314



Introduction

Rabbit anti-DLDH/DLD monoclonal antibody (P09622). For WB, IHC-P, IF/ICC, ELISA. Cross-reacts with Human, Mouse, Rat. Suitable for mitochondrial dehydrogenase complex and histone succinylation studies.

[Learn More](#)

Parameter	Value
Product Name	DLDH/DLD Rabbit mAb
Remarks/Alias	E3; LAD; DLDD; DLDH; GCSL; PHE3; OGDC-E3; DLDH/DLD
Species	Human
GeneID (Human)	1738
GeneID	1738
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 401-500 of human DLDH/DLD (P09622).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09622
Protein Weight	54kDa
Shipping	Ice bag