

Anti-Dldh/Dld Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P09622

Item Number: CM0026331



Introduction

Rabbit polyclonal antibody against human DLDH/DLD, validated for WB, IHC-P, ELISA. Recognizes human and mouse dihydrolipoyl dehydrogenase, a mitochondrial E3 component. 54kDa. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	DLDH/DLD Rabbit pAb
Remarks/Alias	E3; LAD; DLDD; DLDH; GCSL; PHE3; OGDC-E3; DLDH/DLD
Species	Human
GeneID (Human)	1738
GeneID	1738
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 36-300 of human DLDH/DLDH/DLDH/DLD (NP_000099.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09622
Protein Weight	54kDa
Shipping	Ice bag