


Anti-Dlg1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q12959

Item Number: CM0026706



Introduction

Rabbit polyclonal antibody targeting human DLG1 (SAP97) for Western blot, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for research into cell polarity, synaptic biology, and cardiac function.

[Learn More](#)

Parameter	Specification
Product Name	DLG1 Rabbit pAb
Remarks/Alias	hdlg; DLGH1; SAP97; SAP-97; DLG1
Species	Human
Gene ID (Human)	1739
Gene ID	1739
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human DLG1 (NP_004078.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12959
Protein Weight	100kDa
Shipping	Ice bag