


Anti-Dlg3 Rabbit Polyclonal Antibody For Wb - Q92796

Item Number: CM0019901



Introduction

Rabbit polyclonal antibody against human DLG3 (SAP102), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for synaptic plasticity and learning research. SWISS: Q92796.

[Learn More](#)

Parameter	Value
Product Name	DLG3 Rabbit pAb
Remarks/Alias	MRX; XLMR; MRX90; NEDLG; SAP102; XLID90; PPP1R82; DLG3
Species	Human
GeneID (Human)	1741
GeneID	1741
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 464-817 of human DLG3 (NP_066943.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92796
Protein Weight	90kDa
Shipping	Ice bag