

Anti-Dlgap2 Polyclonal Antibody For Wb, Elisa - Q9P1A6

Item Number: CM0018232

Introduction

DLGAP2 Rabbit Polyclonal Antibody validated for WB and ELISA. Targets human DLGAP2, a postsynaptic scaffolding protein involved in synaptic organization. Ideal for neuroscience and signaling research.

[Learn More](#)

Parameter	Value
Product Name	DLGAP2 Rabbit pAb
Remarks/Alias	DAP2; SAPAP2; C8orf68; ERICH1-AS1; DLGAP2
Species	Human
GeneID (Human)	9228
GeneID	9228
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 550-800 of human DLGAP2 (NP_004736.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9P1A6
Protein Weight	118kDa
Shipping	Ice bag