

Anti-Dlk1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P80370

Item Number: CM0026096



Introduction

DLK1 Rabbit pAb is a polyclonal antibody raised against recombinant human DLK1 (aa 24-303). Validated for WB, IHC-P, IF/ICC, ELISA, with cross-reactivity to mouse and rat. Suitable for neuroendocrine and developmental research.

[Learn More](#)

Parameter	Value
Product Name	DLK1 Rabbit pAb
Remarks/Alias	DLK; FA1; ZOG; pG2; DLK-1; PREF1; Delta1; Pref-1; DLK1
Species	Human
GeneID	8788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-303 of human DLK1 (P80370)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P80370
Protein Weight	41kDa
Shipping	Ice bag