

Anti-Dlk1 Rabbit Polyclonal Antibody For Wb, Elisa - P80370

Item Number: CM0027224



Introduction

DLK1 Rabbit Polyclonal Antibody (P80370) validated for WB and ELISA with cross-reactivity to human and mouse. Suitable for studies on neuroendocrine differentiation. Source: Rabbit.

[Learn More](#)

Parameter	Value
Product Name	DLK1 Rabbit pAb
Remarks/Alias	DLK; FA1; ZOG; pG2; DLK-1; PREF1; Delta1; Pref-1; DLK1
Species	Human
Gene ID (Human)	8788
Gene ID	8788
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P80370
Protein Weight	41kDa
Shipping	Ice bag