

Anti-DLL1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O00548

Item Number: CM0017263



Introduction

Rabbit polyclonal antibody targeting human DLL1, a Notch ligand crucial for embryonic development and stem cell maintenance. Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for studying Notch signaling, cell differentiation, and tissue patterning.

[Learn More](#)

Parameter	Value
Product Name	DLL1 Rabbit pAb
Remarks/Alias	DL1; Delta; DELTA1; NEDBAS; DLL1
Species	Human
GeneID (Human)	28514
GeneID	28514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-220 of human DLL1 (NP_005609.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00548
Protein Weight	78kDa
Shipping	Ice bag