

Anti-DLL4 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nr61

Item Number: CM0013376



Introduction

High-quality rabbit polyclonal antibody targeting DLL4, a Notch ligand involved in angiogenesis and development. Detects human, mouse, and rat DLL4 in WB and ELISA. Ideal for vascular biology and neurogenesis research.

[Learn More](#)

Parameter	Value
Product Name	DLL4 Rabbit pAb
Remarks/Alias	AOS6; delta4; hdelta2; DLL4
Species	Human
GeneID (Human)	54567
GeneID	54567
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 317-530 of human DLL4 (NP_061947.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NR61
Protein Weight	75kDa
Shipping	Ice bag