

Anti-DLL4 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nr61

Item Number: CM0019962



Introduction

Rabbit polyclonal antibody against human DLL4, validated for Western blotting and ELISA. Ideal for studying Notch signaling, angiogenesis, and neurogenesis. Cross-reacts with human DLL4.

[Learn More](#)

Parameter	Value
Product Name	DLL4 Rabbit pAb
Remarks/Alias	DLL4, AOS6, hdelta2, delta-like protein 4
Species	Human
GeneID	54567
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-524 of human DLL4 (NP_061947.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NR61
Protein Weight	74kDa
Shipping	Ice bag