


Anti-DLST Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P36957

Item Number: CM0017037



Introduction

Rabbit anti-DLST polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat DLST, a key enzyme in the 2-oxoglutarate dehydrogenase complex. Suitable for mitochondrial and nuclear succinylation studies.

[Learn More](#)

Parameter	Value
Product Name	DLST Rabbit pAb
Remarks/Alias	DLTS; KGD2; PGL7; DLST
Species	Human
GeneID (Human)	1743
GeneID	1743
Immunogen	synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 150-250 of human DLST (NP_001924.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36957
Protein Weight	49kDa
Shipping	Ice bag